

WEB337 - *in vivo* biotransformation

using SOMT-2

Benjamin Weigel

February 14, 2015

Contents

1	Introduction	1
2	Experimental	1
2.1	seed culture	1
2.2	main culture	1
2.3	OD ⁶⁰⁰ measurements	2
2.4	HPLC	2
3	Results	2
3.1	OD ⁶⁰⁰	2
4	Appendix	2

List of Figures

1	OD ⁶⁰⁰ of samples	5
---	--	---

List of Tables

1 Introduction

Test whether different substrates available in-lab are converted by SOMT-2 *in vivo*. Use SOMT seed culture to inoculate main cultures. Add substrates after 4 hours of incubation at 30 °C. 16 substrates means 16 flasks. Take two samples from each flask at 0, 10, 20 and 30 hours. $16 \times 4 \times 2 = 128$ samples.

2 Experimental

2.1 seed culture

- 1) ~10 mL pre-culture in LB supplemented with proper AB (100 ug/mL kanamycin)
- 2) grow over night at 30 °C and 220 rpm

2.2 main culture

- 1) pellet cells (5 min @ 5000×g, 4 °C) and wash with 15 mL PBS
- 2) resuspend pellet in the 3 mL of PBS
!! measure and record OD⁶⁰⁰
- 3) inoculate 200 mL of autoinduction medium (+ 100 ug/mL kan) to an OD⁶⁰⁰ = 0.1
- 4) aliquot 10 mL into new flasks for each sample (17 flasks) (*use 100 mL flasks*)
- 5) add 0.1 mM of flavonoid (see 4) from 10 mM stock in MeOH or DMSO to the cultures at 4 hours after inoculation (OD⁶⁰⁰ ~ 0.8)
- 6) take a 600 µl sample at 10, 20 and 30 hours after inoculation and divide as follows:
(**on ice!**)
 - a) measure OD⁶⁰⁰ (~100 µl)
 - b) 500 µl for HPLC (see 2.4)

2.3 OD⁶⁰⁰ measurements

- measure OD⁶⁰⁰ in MTP (all samples, 100 µl of sample) **!pathlength differs from cuvette!**
- measure OD⁶⁰⁰ of random samples in cuvette as reference

2.4 HPLC

1. extract 500 µl of culture **twice** with 500 µl ethyl acetate + 1% formic acid
2. vortex for 30 s to extract, centrifuge for 10 min @ 10.000×g, 4 °C to separate phases
3. pool organic phases and evaporate in SpeedVac (45 °C)
4. solve remainder in 200 µl MeOH
5. analyze via HPLC

Column:

Injection volume: 10 µl

Solvent A: H₂O + 0.2% formic acid

Solvent B: MeCN + 0.2% formic acid

Program: 5% B (hold 4 min) → 21 min ramp → 100% B (hold 5 min)















Figure 1: OD⁶⁰⁰ of samples.



3 Results

3.1 OD⁶⁰⁰

4 Appendix

No.	substrate	moiety	structure
1	alizarin	anthrachinone	
2	purpurin	anthrachinone	
3	apigenin	flavone	
4	chrysin	flavone	

No.	substrate	moiety	structure
5	genistein	isoflavone	
6	galangin	flavonol	
7	kaempferol	flavonol	
8	quercetin	flavonol	
8	reosmin	homo-dihydro cinnamic keton	
9	naringenin	flavanon	
10	eriodictyol	flavanon	
11	homoeriodictyol	flavanon	
12	hesperetin	flavanon	
13	phloretin	chalcon	
14	resveratrol	stilbene	

No.	substrate	moiety	structure
15	<i>p</i> -coumaric acid	cinnamic acid	
16	caffeic acid	cinnamic acid	
17	none (blank)		